

Product datasheet for **HP206872**

Rad6 (UBE2A) Human qPCR Primer Pair (NM_003336)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7319
Forward Sequence:	CTGGACATACTTCAGAACCGTTG
Reverse Sequence:	GGCTGTTTGCTGGACTATTGGG
ACCN:	<u>NM_003336</u> , <u>NM_003336.1</u> , <u>NM_003336.2</u> , <u>NM_003336.3</u> , <u>BC010175</u> , <u>BC010175.2</u> , <u>BC042021</u> , <u>BI093289</u> , <u>NM_003336.4</u>
UniProt ID:	<u>P49459</u>
Synonyms:	HHR6A; MRXS30; MRXSN; RAD6A; UBC2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)