

Product datasheet for **HP206873**

UBE2B Human qPCR Primer Pair (NM_003337)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7320
Forward Sequence:	TCCTTCAGAATCGATGGAGTCCA
Reverse Sequence:	GCTGGACTGTTAGGATTCGGTTC
ACCN:	<u>BC001694</u> , <u>NM_003337</u> , <u>NM_003337.1</u> , <u>NM_003337.2</u> , <u>NM_003337.3</u> , <u>BC005979</u> , <u>BC008404</u> , <u>BC008470</u> , <u>BT007071</u> , <u>BU954901</u> , <u>NM_003337.4</u>
UniProt ID:	<u>P63146</u>
Synonyms:	E2-17kDa; HHR6B; HR6B; RAD6B; UBC2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)