

Product datasheet for **HP206875**

UBE2D2 Human qPCR Primer Pair (NM_003339)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7322
Forward Sequence:	CTACGATCACAGTGGTCTCCAG
Reverse Sequence:	CGAGCAATCTCAGGCACTAAAGG
ACCN:	<u>NM_003339</u> , <u>NM_003339.1</u> , <u>NM_003339.2</u> , <u>BC033349</u> , <u>BC033349.1</u> , <u>BC019052</u> , <u>BC048426</u> , <u>BQ186099</u> , <u>NM_003339.3</u>
UniProt ID:	<u>P62837</u>
Synonyms:	E2(17)KB2; PUBC1; UBC4; UBC4/5; UBCH4; UBCH5B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)