

Product datasheet for **HP206885**

MMS2 (UBE2V2) Human qPCR Primer Pair (NM_003350)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7336
Forward Sequence:	AAGGAGTAGGCGACGGTACAGT
Reverse Sequence:	AGGTCCACATTCTACTTTCAGGC
ACCN:	<u>NM_003350</u> , <u>NM_003350.2</u> , <u>BC016710</u> , <u>BC016710.1</u> , <u>BC007051</u> , <u>BC016332</u> , <u>BC028673</u> , <u>BC062418</u> , <u>BT006744</u> , <u>NM_003350.3</u>
UniProt ID:	<u>Q15819</u>
Synonyms:	DDVit-1; DDVIT1; EDAF-1; EDPF-1; EDPF1; MMS2; UEV-2; UEV2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)