

Product datasheet for **HP206900**

VAV2 Human qPCR Primer Pair (NM_003371)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 7410 |
| Forward Sequence: | CTGCTGTTCCACAAGATGACCG |
| Reverse Sequence: | AGCCCTGCTTTCCTTGAAGTG |
| ACCN: | NM_003371 , NM_003371.1 , NM_003371.2 , NM_003371.3 , BC132965 , BC033187 , BC043615 , BC128448 , BC132967 , BQ050947 , BQ071939 , BX640754 , NM_003371.4 |
| UniProt ID: | P52735 |
| Synonyms: | VAV-2 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)