

## Product datasheet for **HP206913**

### CYLN2 (CLIP2) Human qPCR Primer Pair (NM\_003388)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7461
Forward Sequence:	ACATCTGCGAGGTGGAGAAGGA
Reverse Sequence:	AGGTCCACCTTCTCTTCCGCA
ACCN:	<u><a href="#">NM_003388</a></u> , <u><a href="#">NM_003388.2</a></u> , <u><a href="#">NM_003388.3</a></u> , <u><a href="#">BC006259</a></u> , <u><a href="#">BC140369</a></u> , <u><a href="#">BC148702</a></u> , <u><a href="#">BC150332</a></u> , <u><a href="#">NM_003388.5</a></u>
UniProt ID:	<u><a href="#">Q9UDT6</a></u>
Synonyms:	CLIP; CLIP-115; CYLN2; WBSCR3; WBSCR4; WSCR3; WSCR4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)