

Product datasheet for **HP206923**

CRM1 (XPO1) Human qPCR Primer Pair (NM_003400)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7514
Forward Sequence:	CTACATCTGCCTCTCCGTTGCT
Reverse Sequence:	CCAATACTTCCTCTGGTTTAGCC
ACCN:	NM_003400 , NM_003400.1 , NM_003400.2 , NM_003400.3 , BC032847 , BC032847.2 , BX492941 , BX643790 , BX647758 , BX954288 , NM_003400.4
UniProt ID:	O14980
Synonyms:	CRM-1; CRM1; emb; exp1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)