

## Product datasheet for **HP206943**

### Myeloid zinc finger 1 (MZF1) Human qPCR Primer Pair (NM\_003422)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7593
Forward Sequence:	GCAGGTGAAAGAGGAGTCAGAG
Reverse Sequence:	GTCCAGCACGGTCCCACATCT
ACCN:	<a href="#">BC007777</a> , <a href="#">NM_003422</a> , <a href="#">NM_003422.1</a> , <a href="#">NM_003422.2</a> , <a href="#">BC007777.2</a> , <a href="#">BC053316</a> , <a href="#">BQ446552</a> , <a href="#">BU620994</a>
UniProt ID:	<a href="#">P28698</a>
Synonyms:	MZF-1; MZF1B; ZFP98; ZNF42; ZSCAN6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)