

Product datasheet for **HP206947**

ZNF76 Human qPCR Primer Pair (NM_003427)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7629
Forward Sequence:	GTCATACAGGTGACCGTCCAT
Reverse Sequence:	GGTTTCTACCAAGTGTGGGTAC
ACCN:	NM_003427 , NM_003427.1 , NM_003427.2 , NM_003427.3 , NM_003427.4 , BC002549 , BC002549.2 , BC000133 , BE962441 , BT007316
UniProt ID:	P36508
Synonyms:	D6S229E; Zfp523; ZNF523
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)