

Product datasheet for **HP206970**

ZFP200 (ZNF200) Human qPCR Primer Pair (NM_003454)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7752
Forward Sequence:	TAGTGCTGGAGCACTTCCTCAC
Reverse Sequence:	TTCACCAAGGGACGAGGATGCA
ACCN:	NM_003454 , NM_003454.1 , NM_003454.2 , NM_003454.3 , BC032575 , BC032575.1 , BC054005 , BC012909 , BC025293 , BI459882 , BI822625 , BU177850 , NM_003454.4
UniProt ID:	P98182
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)