

## Product datasheet for **HP206976**

### PTP4A1 Human qPCR Primer Pair (NM\_003463)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7803
Forward Sequence:	AGGTGGAATGAAATACGAAGATGC
Reverse Sequence:	GGAATCTTTGAAACGCAGCCGC
ACCN:	<a href="#">NM_003463</a> , <a href="#">NM_003463.1</a> , <a href="#">NM_003463.2</a> , <a href="#">NM_003463.3</a> , <a href="#">NM_003463.4</a> , <a href="#">BC023975</a> , <a href="#">BC023975.2</a> , <a href="#">BC045571</a> , <a href="#">BF436021</a> , <a href="#">BG717869</a> , <a href="#">BI222469</a> , <a href="#">BP264358</a> , <a href="#">BU632706</a>
UniProt ID:	<a href="#">Q93096</a>
Synonyms:	HH72; PRL-1; PRL1; PTP(CAAX1); PTPCAAX1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)