

Product datasheet for **HP206980**

Chromogranin C (SCG2) Human qPCR Primer Pair (NM_003469)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7857
Forward Sequence:	GAGAAGCCGAATGGATCAGTGG
Reverse Sequence:	TCTGGATGGTCTAAGTCAGCCTC
ACCN:	NM_003469 , NM_003469.1 , NM_003469.2 , NM_003469.3 , NM_003469.4 , BC022509 , BC022509.1 , NM_003469.5
UniProt ID:	P13521
Synonyms:	CHGC; EM66; SgII; SN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)