

## Product datasheet for **HP206983**

### DEK Human qPCR Primer Pair (NM\_003472)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7913
Forward Sequence:	TGGGTCAGTTCAGTGGCTTTCC
Reverse Sequence:	CTCTCAAATCAAGAACCTCACAG
ACCN:	<a href="#">NM_003472</a> , <a href="#">NM_003472.1</a> , <a href="#">NM_003472.2</a> , <a href="#">NM_003472.3</a> , <a href="#">BC035259</a> , <a href="#">BC035259.1</a> , <a href="#">BC012442</a> , <a href="#">BQ574048</a> , <a href="#">BX641063</a> , <a href="#">NM_003472.4</a>
UniProt ID:	<a href="#">P35659</a>
Synonyms:	D6S231E
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)