

Product datasheet for **HP206997**

RIP140 (NRIP1) Human qPCR Primer Pair (NM_003489)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 8204 |
| Forward Sequence: | AGCAAAGACCCACCAGGAGAGA |
| Reverse Sequence: | GTCTCTGCTCTTCCACTGACATG |
| ACCN: | <u>NM_003489</u> , <u>NM_003489.1</u> , <u>NM_003489.2</u> , <u>NM_003489.3</u> , <u>BC040361</u> , <u>BE792046</u> , <u>NM_003489.4</u> |
| UniProt ID: | <u>P48552</u> |
| Synonyms: | CAKUT3; RIP140 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)