

Product datasheet for **HP206999**

ARD1A (NAA10) Human qPCR Primer Pair (NM_003491)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8260
Forward Sequence:	TGCTGAGGACGAGAATGGGAAG
Reverse Sequence:	CTGGTCCATCAGTTTCTGAGCC
ACCN:	<u>NM_003491</u> , <u>NM_003491.1</u> , <u>NM_003491.2</u> , <u>NM_003491.3</u> , <u>BC019312</u> , <u>BC019312.1</u> , <u>BC000308</u> , <u>BC063377</u> , <u>BI192945</u> , <u>NM_003491.4</u>
UniProt ID:	<u>P41227</u>
Synonyms:	ARD1; ARD1A; ARD1P; DXS707; hARD1; MCOPS1; NATD; OGDNS; TE2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)