

Product datasheet for **HP207005**

Stannin (SNN) Human qPCR Primer Pair (NM_003498)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8303
Forward Sequence:	AGTCAGAGGACGAGGAGAGCAT
Reverse Sequence:	ATCCTGGCTCAGCCGTGGACT
ACCN:	<u>NM_003498</u> , <u>NM_003498.1</u> , <u>NM_003498.2</u> , <u>NM_003498.3</u> , <u>NM_003498.4</u> , <u>NM_003498.5</u> , <u>BC036443</u> , <u>BC036443.1</u> , <u>BC005837</u> , <u>BC036100</u> , <u>BQ019000</u> , <u>NM_003498.6</u>
UniProt ID:	<u>O75324</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)