

## Product datasheet for **HP207007**

### Axin 1 (AXIN1) Human qPCR Primer Pair (NM\_003502)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8312
Forward Sequence:	GTATGTGCAGGAGGTTATGCCG
Reverse Sequence:	CACCTTCCTCTGCGATCTTGTC
ACCN:	<a href="#">NM_003502</a> , <a href="#">NM_003502.1</a> , <a href="#">NM_003502.2</a> , <a href="#">NM_003502.3</a> , <a href="#">BC044648</a> , <a href="#">BC044648.1</a> , <a href="#">BC017447</a> , <a href="#">BC035872</a> , <a href="#">BX648449</a> , <a href="#">NM_003502.4</a>
UniProt ID:	<a href="#">O15169</a>
Synonyms:	AXIN; PPP1R49
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)