

Product datasheet for **HP207022**

H2BC8 Human qPCR Primer Pair (NM_003518)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8339
Forward Sequence:	CGATACTGGCATCTCATCCAAGG
Reverse Sequence:	ATCTCCCTGGAGGTAATGGTCG
ACCN:	<u>NM_003518</u> , <u>NM_003518.1</u> , <u>NM_003518.2</u> , <u>NM_003518.3</u> , <u>BC082232</u> , <u>NM_003518.4</u>
UniProt ID:	<u>P62807</u>
Synonyms:	dj221C16.8; H2B.1A; H2B/a; H2BC4; H2BC6; H2BC7; H2BC10; H2BFA; HIST1H2BG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)