

## Product datasheet for **HP207025**

### H2BC14 Human qPCR Primer Pair (NM\_003521)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8342
Forward Sequence:	GTACAAGGTGCTGAAGCAGGTC
Reverse Sequence:	ATCTCCCTCGAAGTGATGGTCG
ACCN:	<a href="#">NM_003521</a> , <a href="#">NM_003521.1</a> , <a href="#">NM_003521.2</a> , <a href="#">BC066244</a> , <a href="#">BC067485</a> , <a href="#">BC067486</a> , <a href="#">BC067487</a> , <a href="#">BC067488</a> , <a href="#">BC067489</a> , <a href="#">NM_003521.3</a>
UniProt ID:	<a href="#">Q99879</a>
Synonyms:	dj160A22.3; H2B/e; H2BFE; HIST1H2BM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)