

Product datasheet for **HP207030**

H2BC4 Human qPCR Primer Pair (NM_003526)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 8347 |
| Forward Sequence: | GACACTGGCATCTCTTCCAAGG |
| Reverse Sequence: | AGGTGATGGTCGAGCGCTTGTT |
| ACCN: | NM_003526 , NM_003526.1 , NM_003526.2 , BC009612 , BC105595 , BC106899 , BC160109 , NM_003526.3 |
| UniProt ID: | P62807 |
| Synonyms: | dj221C16.3; H2B.1; H2B/l; H2BC6; H2BC7; H2BC8; H2BC10; H2BFL; HIST1H2BC |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)