

## Product datasheet for **HP207031**

### H2BC17 Human qPCR Primer Pair (NM\_003527)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8348
Forward Sequence:	CAAGGTGCTGAAGCAAGTCCAC
Reverse Sequence:	AGGTGATGGTCGAGCGCTTGTT
ACCN:	<u><a href="#">NM_003527</a></u> , <u><a href="#">NM_003527.1</a></u> , <u><a href="#">NM_003527.2</a></u> , <u><a href="#">NM_003527.3</a></u> , <u><a href="#">NM_003527.4</a></u> , <u><a href="#">BC106720</a></u> , <u><a href="#">BX647290</a></u>
UniProt ID:	<u><a href="#">P23527</a></u>
Synonyms:	dj193B12.2; H2B.2; H2B/n; H2BFN; HIST1H2BO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)