

Product datasheet for **HP207033**

H3FA (HIST1H3A) Human qPCR Primer Pair (NM_003529)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8350
Forward Sequence:	TCCGCCGTTATCAGAAGTCCAC
Reverse Sequence:	GCTCTGGAAACGCAGGTCTGTT
ACCN:	NM_003529 , NM_003529.1 , NM_003529.2 , BC066245 , BC066246 , BC066247 , BC067490 , BC067491 , BC069303 , NM_003529.3
UniProt ID:	P68431
Synonyms:	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3C10; H3C11; H3C12; H3FA; HIST1H3A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)