

Product datasheet for **HP207036**

H3FD (HIST1H3E) Human qPCR Primer Pair (NM_003532)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8353
Forward Sequence:	GGTGCGAGAAATAGCTCAGGAC
Reverse Sequence:	GAATAGCGCACAGTTGGTGTC
ACCN:	<u>NM_003532</u> , <u>NM_003532.1</u> , <u>NM_003532.2</u> , <u>NM_003532.3</u>
UniProt ID:	<u>P68431</u>
Synonyms:	H3.1; H3/d; H3C1; H3C2; H3C3; H3C4; H3C7; H3C8; H3C10; H3C11; H3C12; H3FD; HIST1H3E
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)