

## Product datasheet for **HP207039**

### H3FJ (HIST1H3J) Human qPCR Primer Pair (NM\_003535)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8356
Forward Sequence:	TCCGCAAAGTCCATTTTCAGCG
Reverse Sequence:	GGCGTGAATAGCACAGAGGTTG
ACCN:	<u><a href="#">NM_003535</a></u> , <u><a href="#">NM_003535.1</a></u> , <u><a href="#">NM_003535.2</a></u> , <u><a href="#">BC111921</a></u> , <u><a href="#">NM_003535.3</a></u>
UniProt ID:	<u><a href="#">P68431</a></u>
Synonyms:	H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3C10; H3C11; H3FJ; HIST1H3J
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)