

Product datasheet for **HP207042**

H4C1 Human qPCR Primer Pair (NM_003538)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8359
Forward Sequence:	TGTTGCGTGACAACATCCAGGG
Reverse Sequence:	CTCCGTATAGGTGACAGCGTCA
ACCN:	NM_003538 , NM_003538.1 , NM_003538.2 , NM_003538.3 , BC069654 , BC066248 , BC066249 , BC067495 , NM_003538.4
UniProt ID:	P62805
Synonyms:	H4-16; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C14; H4C15; H4FA; HIST1H4A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)