

Product datasheet for **HP207043**

H4C4 Human qPCR Primer Pair (NM_003539)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8360
Forward Sequence:	AGGAATCACCAAGCCCGCTATC
Reverse Sequence:	TTCCGTGTAGGTGACAGCATCG
ACCN:	<u>NM_003539</u> , <u>NM_003539.1</u> , <u>NM_003539.2</u> , <u>NM_003539.3</u> , <u>BC128104</u> , <u>BC128105</u> , <u>BC143045</u> , <u>BC156842</u> , <u>NM_003539.4</u>
UniProt ID:	<u>P62805</u>
Synonyms:	dj221C16.9; H4-16; H4/b; H4C1; H4C2; H4C3; H4C5; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C14; H4C15; H4FB; HIST1H4D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)