

Product datasheet for **HP207064**

ULK1 Human qPCR Primer Pair (NM_003565)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8408
Forward Sequence:	GCAAGGACTCTTCCTGTGACAC
Reverse Sequence:	CCACTGCACATCAGGCTGTCTG
ACCN:	<u>NM_003565</u> , <u>NM_003565.1</u> , <u>NM_003565.2</u> , <u>BC111603</u> , <u>BI822269</u> , <u>BU729635</u> , <u>NM_003565.4</u>
UniProt ID:	<u>O75385</u>
Synonyms:	ATG1; ATG1A; hATG1; UNC51; Unc51.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)