

Product datasheet for **HP207128**

DGKE Human qPCR Primer Pair (NM_003647)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8526
Forward Sequence:	TGGGTCTTTCCACTGTGCTCAG
Reverse Sequence:	TCCACCTGCATTGGCATCATGG
ACCN:	NM_003647 , NM_003647.1 , NM_003647.2 , BC022297 , BC130629 , BC130631 , BM682264 , NM_003647.3
UniProt ID:	P52429
Synonyms:	AHUS7; DAGK5; DAGK6; DGK; NPHS7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)