

Product datasheet for **HP207138**

BARX2 Human qPCR Primer Pair (NM_003658)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8538
Forward Sequence:	ACTTGGCTCAGTCTCTGGGACT
Reverse Sequence:	GGTTTTGTGGGTGCTTCCTGTC
ACCN:	NM_003658 , NM_003658.1 , NM_003658.2 , NM_003658.3 , NM_003658.4 , BC069378 , BC069378.1 , BC111432 , BC111572 , NM_003658.5
UniProt ID:	Q9UMQ3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)