

## Product datasheet for **HP207140**

### PPFIA3 Human qPCR Primer Pair (NM\_003660)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8541
Forward Sequence:	CGATGACCACAATAAGCGGCTG
Reverse Sequence:	CTTCATGTTGGCTATCTCCTCGC
ACCN:	<u><a href="#">NM_003660</a></u> , <u><a href="#">NM_003660.1</a></u> , <u><a href="#">NM_003660.2</a></u> , <u><a href="#">NM_003660.3</a></u> , <u><a href="#">BC101518</a></u> , <u><a href="#">BC101518.1</a></u> , <u><a href="#">BC021255</a></u> , <u><a href="#">BC101520</a></u> , <u><a href="#">BC143470</a></u> , <u><a href="#">BU632015</a></u> , <u><a href="#">BX111585</a></u> , <u><a href="#">NM_003660.4</a></u>
UniProt ID:	<u><a href="#">O75145</a></u>
Synonyms:	LPNA3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)