

Product datasheet for **HP207147**

INE1 Human qPCR Primer Pair (NM_003669)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8552
Forward Sequence:	CTTCTGCTGTGGAAGAGCCTGC
Reverse Sequence:	TGGCACACAGACAAGCAGGATC
ACCN:	<u>NM_003669</u> , <u>NM_003669.1</u> , <u>BC069671</u> , <u>BC069696</u> , <u>BC069772</u> , <u>BC112149</u> , <u>BC112175</u>
Synonyms:	MGC97326, MGC97335, MGC97386
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)