

Product datasheet for **HP207151**

PRP18 homolog (PRPF18) Human qPCR Primer Pair (NM_003675)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8559
Forward Sequence:	GGTATCCATGCCAGAACTGGCA
Reverse Sequence:	GGAAAGTGTCTTGCCAAATGGTC
ACCN:	<u>BC002572, NM_003675, NM_003675.1, NM_003675.2, NM_003675.3, BC000794, NM_003675.4</u>
UniProt ID:	<u>Q99633</u>
Synonyms:	hPrp18; PRP18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)