

## Product datasheet for **HP207163**

### CASK Human qPCR Primer Pair (NM\_003688)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8573
Forward Sequence:	CAGGAATCCTCCAAGCGATGCA
Reverse Sequence:	CTTCATGTGCCACTACGTCGTG
ACCN:	<u><a href="#">NM_003688</a></u> , <u><a href="#">NM_003688.1</a></u> , <u><a href="#">NM_003688.2</a></u> , <u><a href="#">NM_003688.3</a></u> , <u><a href="#">BC117311</a></u> , <u><a href="#">BC029936</a></u> , <u><a href="#">BC035990</a></u> , <u><a href="#">BC143454</a></u> , <u><a href="#">BC143456</a></u>
UniProt ID:	<u><a href="#">O14936</a></u>
Synonyms:	CAGH39; CAMGUK; CMG; FGS4; hCASK; LIN2; MICPCH; MRXSNA; TNRC8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)