

Product datasheet for **HP207175**

FAM193A Human qPCR Primer Pair (NM_003704)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8603
Forward Sequence:	CCAACCAGAAGGAGGGCAAGTA
Reverse Sequence:	CTCTTGGTCAGCCGTAAGCGAA
ACCN:	<u>NM_003704</u> , <u>NM_003704.1</u> , <u>NM_003704.2</u> , <u>NM_003704.3</u> , <u>BC136646</u> , <u>BC144337</u> , <u>BC144339</u> , <u>BC144340</u> , <u>BF924432</u> , <u>BM468825</u> , <u>BM788933</u>
UniProt ID:	<u>P78312</u>
Synonyms:	C4orf8; RES4-22
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)