

Product datasheet for **HP207198**

OASL Human qPCR Primer Pair (NM_003733)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8638
Forward Sequence:	GTGCCTGAAACAGGACTGTTGC
Reverse Sequence:	CCTCTGCTCCACTGTCAAGTGG
ACCN:	<u>NM_003733</u> , <u>NM_003733.1</u> , <u>NM_003733.2</u> , <u>NM_003733.3</u> , <u>BC117408</u> , <u>BC117410</u> , <u>BG339739</u> , <u>BG744291</u> , <u>BM909951</u> , <u>BM918192</u>
UniProt ID:	<u>Q15646</u>
Synonyms:	OASL1; OASLd; p59 OASL; p59-OASL; p59OASL; TRIP-14; TRIP14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)