

Product datasheet for **HP207208**

NCOA1 Human qPCR Primer Pair (NM_003743)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8648
Forward Sequence:	ATCGGAGATCCTGCCAGCTTCA
Reverse Sequence:	ATGCCTGTTCTGGTTGTCCA
ACCN:	<u>NM_003743</u> , <u>NM_003743.1</u> , <u>NM_003743.2</u> , <u>NM_003743.3</u> , <u>NM_003743.4</u> , <u>BC111533</u> , <u>BC111534</u> , <u>BC035976</u> , <u>BE675986</u> , <u>BF975824</u> , <u>NM_003743.5</u>
UniProt ID:	<u>Q15788</u>
Synonyms:	bHLHe42; bHLHe74; F-SRC-1; KAT13A; RIP160; SRC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)