

Product datasheet for **HP207216**

EIF3S4 (EIF3G) Human qPCR Primer Pair (NM_003755)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8666
Forward Sequence:	CGCATCTACCTGGCTAAGGACA
Reverse Sequence:	AGGATGAGGTGGTCGTAGCCAA
ACCN:	NM_003755 , NM_003755.1 , NM_003755.2 , NM_003755.3 , NM_003755.4 , BC008469 , BC008469.1 , BC000733 , BE017996 , BF701069 , BT006889 , NM_003755.5
UniProt ID:	O75821
Synonyms:	eIF3-delta; EIF3-P42; eIF3-p44; EIF3S4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)