

Product datasheet for **HP207219**

EIF3S1 (EIF3J) Human qPCR Primer Pair (NM_003758)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8669
Forward Sequence:	TGGAAGTCTTAGTTCGAGATGTGT
Reverse Sequence:	CCTCCAGGAACACACCTTTCT
ACCN:	<u>NM_003758</u> , <u>NM_003758.1</u> , <u>NM_003758.2</u> , <u>NM_003758.3</u> , <u>BC002719</u> , <u>BC002719.2</u> , <u>BU621047</u> , <u>NM_003758.4</u>
UniProt ID:	<u>O75822</u>
Synonyms:	eIF3-alpha; eIF3-p35; EIF3S1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)