

Product datasheet for **HP207223**

STX11 Human qPCR Primer Pair (NM_003764)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8676
Forward Sequence:	GAGATGAAGCAGCGGACAACCT
Reverse Sequence:	CCAGCAAGTTCTCGGAAAACACG
ACCN:	NM_003764 , NM_003764.1 , NM_003764.2 , NM_003764.3 , BC033519 , BC033519.1 , BF941242 , BM150302 , BQ030324 , NM_003764.4
UniProt ID:	O75558
Synonyms:	FHL4; HLH4; HPLH4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)