

Product datasheet for **HP207234**

MRPL40 Human qPCR Primer Pair (NM_003776)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64976
Forward Sequence:	GACCAAGAAGCAAAGGAGCGCT
Reverse Sequence:	CCTCTCAGTCTCCTCAAAGGTG
ACCN:	<u>NM_003776</u> , <u>NM_003776.1</u> , <u>NM_003776.2</u> , <u>NM_003776.3</u> , <u>BC009707</u> , <u>BM974354</u> , <u>BP309678</u> , <u>BX382916</u>
UniProt ID:	<u>Q9NQ50</u>
Synonyms:	L40mt; MRP-L22; MRP-L40; MRPL22; NLVCF; URIM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)