

Product datasheet for **HP207248**

SNX3 Human qPCR Primer Pair (NM_003795)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8724
Forward Sequence:	ACAGTGACTTTGAATGGCTGCGA
Reverse Sequence:	CTCCTCTAAAAGGAAGCTGACGC
ACCN:	<u>NM_003795</u> , <u>NM_003795.1</u> , <u>NM_003795.2</u> , <u>NM_003795.3</u> , <u>NM_003795.4</u> , <u>NM_003795.5</u> , <u>BC015179</u> , <u>BC015179.1</u> , <u>BC008444</u> , <u>BC014580</u> , <u>BC016863</u> , <u>BG497532</u> , <u>BQ432490</u> , <u>BT007114</u> , <u>NM_003795.6</u>
UniProt ID:	<u>O60493</u>
Synonyms:	Grd19; MCOPS8; SDP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)