

Product datasheet for **HP207252**

HCM (RNMT) Human qPCR Primer Pair (NM_003799)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8731
Forward Sequence:	CTGTTTTGGACCTGGGATGTGG
Reverse Sequence:	TCCTCATACCGCTGCTGACACT
ACCN:	<u>NM_003799</u> , <u>NM_003799.1</u> , <u>NM_003799.2</u> , <u>BC036798</u> , <u>BC036798.1</u> , <u>BC005290</u> , <u>BC017816</u> , <u>BX647977</u> , <u>NM_003799.3</u>
UniProt ID:	<u>O43148</u>
Synonyms:	cm1p; CMT1; CMT1c; hCMT1; hMet; MET; Met; RG7MT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)