

Product datasheet for **HP207253**

HCE (RNGTT) Human qPCR Primer Pair (NM_003800)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8732
Forward Sequence:	GGCTTCAATCGCACTGGTTCC
Reverse Sequence:	CCTGTAGATTCTGGTGGTCTG
ACCN:	NM_003800 , NM_003800.1 , NM_003800.2 , NM_003800.3 , NM_003800.4 , BC019954 , BC019954.1 , BP256840 , BP396964 , BX537450 , NM_003800.5
UniProt ID:	O60942
Synonyms:	CAP1A; hCAP; HCE; HCE1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)