

Product datasheet for **HP207254**

GPAA1 Human qPCR Primer Pair (NM_003801)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8733
Forward Sequence:	CCTGCTCAATCTCTCCAGACC
Reverse Sequence:	CAGCGGTCCATCCAATGATGTC
ACCN:	<u>BC003171</u> , <u>BC006383</u> , <u>NM_003801</u> , <u>NM_003801.1</u> , <u>NM_003801.2</u> , <u>NM_003801.3</u> , <u>BC003171.1</u> , <u>BC006383.2</u> , <u>BC004129</u> , <u>BC004129.2</u> , <u>NM_003801.4</u>
UniProt ID:	<u>O43292</u>
Synonyms:	GAA1; GPIBD15; hGAA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)