

## Product datasheet for **HP207270**

### CDS2 Human qPCR Primer Pair (NM\_003818)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8760
Forward Sequence:	CTCCACCTTTGCCTCGCTCATT
Reverse Sequence:	CTGGCAGTCAAAGCGATCCATG
ACCN:	<a href="#">NM_003818</a> , <a href="#">NM_003818.1</a> , <a href="#">NM_003818.2</a> , <a href="#">NM_003818.3</a> , <a href="#">BC025751</a> , <a href="#">BC025751.1</a> , <a href="#">BC015432</a> , <a href="#">BC025346</a> , <a href="#">BC044234</a> , <a href="#">BP353230</a> , <a href="#">NM_003818.4</a>
UniProt ID:	<a href="#">O95674</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)