

## Product datasheet for **HP207277**

### SNAP23 Human qPCR Primer Pair (NM\_003825)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8773
Forward Sequence:	GGTGACAAATGGTCAGCTTCAGC
Reverse Sequence:	CCAGGATACTGCCCACTTGAGT
ACCN:	<a href="#">NM_003825</a> , <a href="#">NM_003825.1</a> , <a href="#">NM_003825.2</a> , <a href="#">NM_003825.3</a> , <a href="#">BC003686</a> , <a href="#">BC003686.1</a> , <a href="#">BC000148</a> , <a href="#">BC022890</a> , <a href="#">BI918294</a> , <a href="#">BM782600</a> , <a href="#">BT006916</a> , <a href="#">NM_003825.4</a>
UniProt ID:	<a href="#">O00161</a>
Synonyms:	HsT17016; SNAP-23; SNAP23A; SNAP23B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)