

Product datasheet for **HP207290**

DcR1 (TNFRSF10C) Human qPCR Primer Pair (NM_003841)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8794
Forward Sequence:	GGTGTGGATTACACCAACGCTTC
Reverse Sequence:	CTGACACACTGTGTCTCTGGTC
ACCN:	<u>NM_003841</u> , <u>NM_003841.1</u> , <u>NM_003841.2</u> , <u>NM_003841.3</u> , <u>NM_003841.4</u> , <u>BC017852</u> , <u>BC125041</u> , <u>BC125042</u> , <u>BG545631</u>
UniProt ID:	<u>O14798</u>
Synonyms:	CD263; DCR1; DCR1-TNFR; LIT; TRAIL-R3; TRAILR3; TRID
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)