

## Product datasheet for **HP207297**

### SUCLG2 Human qPCR Primer Pair (NM\_003848)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 8801   |
| Forward Sequence: | ATGGTGGGAAGCCAGCCAACTT   |
| Reverse Sequence: | TGGCACAGTTGACGATACCACC   |
| ACCN:             | <u><a href="#">NM_003848</a></u> , <u><a href="#">NM_003848.1</a></u> , <u><a href="#">NM_003848.2</a></u> , <u><a href="#">NM_003848.3</a></u> , <u><a href="#">BC047024</a></u> , <u><a href="#">BC007716</a></u> , <u><a href="#">BC019868</a></u> ,<br><u><a href="#">BC035149</a></u> , <u><a href="#">BC068602</a></u>                 |
| UniProt ID:       | <u><a href="#">Q96I99</a></u>  |
| Synonyms:         | G-SCS; GBETA; GTPSCS   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)