

Product datasheet for **HP207326**

CCN5 Human qPCR Primer Pair (NM_003881)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8839
Forward Sequence:	CCTGCGACCAACTCCACGTCT
Reverse Sequence:	TTCACCTCACAGTGCTGTCGT
ACCN:	NM_003881 , NM_003881.1 , NM_003881.2 , NM_003881.3 , BC017782 , BC017782.2 , BC058074 , BC064379 , BI823598 , BU623413
UniProt ID:	O76076
Synonyms:	CT58; CTGF-L; WISP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)