

## Product datasheet for **HP207339**

### IEX1 (IER3) Human qPCR Primer Pair (NM\_003897)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8870
Forward Sequence:	GCAGCCGCAGGGTCTCTACC
Reverse Sequence:	CTCTTCAGCCATCAGGATCTGG
ACCN:	<u><a href="#">NM_003897</a></u> , <u><a href="#">NM_003897.1</a></u> , <u><a href="#">NM_003897.2</a></u> , <u><a href="#">NM_003897.3</a></u> , <u><a href="#">BC005080</a></u> , <u><a href="#">BC000844</a></u> , <u><a href="#">BM915952</a></u> , <u><a href="#">BT006703</a></u>
UniProt ID:	<u><a href="#">P46695</a></u>
Synonyms:	DIF-2; DIF2; GLY96; IEX-1; IEX-1L; IEX1; PRG1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)